



Guth Laboratories^{INC}

Model G-590 Model 590 Simulator System



CLICK ON MODULE NAME IN LEFT HAND COLUMN OR ANY SYSTEM COMPONENT HIGHLIGHTED IN RED FOR DETAILED MODULE INFORMATION.

MODEL 590 COMPLETE PORTABLE SYSTEM WITH AIR PUMP, BATTERY BACK, 12 VOLT DC PLUG, AC CHARGER, RS232 INTERFACE CABLE

MODEL 590 WITH 12 VOLT DC BATTERY PACK MODULE RECHARGEABLE BATTERY BATTERY PACK WITH AN AVERAGE TIME PER CHARGE OF 2 - 3 HOURS. CAPABLE OF PERFORMING 50 TESTS WHEN CONNECTED TO BENCH TOP INSTRUMENTS WITH RECIRCULATION SUPPLIED WITH 110V AC CHARGER.

MODEL 590 WITH PORTABLE SIMULATOR AIR PUMP MODULE PORTABLE PUMP DESIGNED TO BE USED WITH IGNITION INTERLOCK SYSTEMS AND OR PBT'S FOR TESTING AND CALIBRATION PURPOSES. WHEN THE PUMP MODULE IS CONNECTED TO THE MODEL 590 IT WILL DELIVER A SAMPLE OF ALCOHOL VAPOR AT 10 LPM. THE OPERATOR CONTROLS FLOW TIME BY USING THE RUN BUTTON ON THE PUMP BOX.

MODEL 590 WITH AC POWER ADAPTER AC POWER SUPPLY FOR CONTINUOUS POWER FOR THOSE WHO WISH TO USE THE MODEL 590 AS A "BENCH TOP" SIMULATOR.