



# Guth Laboratories<sup>INC</sup>

## Model 590-PS 590 AC POWER ADAPTER



The 590-PS AC power supply for continuous power for those who wish to use the Model 590 as a "bench top" simulator.

### **SPECIFICATIONS**

Construction: All Metal, Powder Coated

Weight: < 1.5 lbs

Size: 7" L x 3 3/4" W x 1.75" H

Output: 13 VDC @ 5 AMPS

Input: 90-264 VAC, 43-63 Hz.

Protection: Fused Input & Output

Warranty: 1 Year Parts & Labor